

---

**neuroI**

***Release 0.0.5***

**Feb 10, 2022**




---

## Contents:

---

|          |  |           |
|----------|--|-----------|
| <b>1</b> | <b>neuroI package</b>                  | <b>3</b>  |
| 1.1      | neuroI.models subpackage . . . . .     | 3         |
| 1.2      | neuroI.BCI module . . . . .            | 9         |
| 1.3      | neuroI.BCI_tools module . . . . .      | 13        |
| 1.4      | neuroI.plot module . . . . .           | 14        |
| 1.5      | neuroI.connect_device module . . . . . | 15        |
| 1.6      | neuroI.streams . . . . .               | 15        |
| <b>2</b> | <b>Installation</b>                    | <b>17</b> |
| <b>3</b> | <b>Indices</b>                         | <b>19</b> |
|          | <b>Index</b>                           | <b>21</b> |





Welcome! `neuro1` is a python package for implementing Brain-Computer Interfaces in a modular manner. With the help of tools in the package, you will be able define the behavior of your intended BCI and easily implement it. A `neuro1` BCI is defined by a number of components:

- A `classifier` which decodes brain data into some kind of 'brain-state'
- An `action` which provides feedback depending on the decoded 'brain-state'
- An optional `calibrator` which runs at startup and modifies the operation of the BCI
- An optional `transformer` which transforms the current `buffer` of data into the form expected by the `classifier`

The `neuro1` BCI manages an incoming stream of brain data and uses the above user-defined functions to run a brain-computer interface.

The package includes generic utility functions to aid in creating the `classifier`'s, `transformer`'s, and `calibrator`'s for common BCI use-cases. It also comes prepackaged with a growing list of trained machine learning models for common BCI classification tasks.



## 1.1 neurol.models subpackage

---

### 1.1.1 neurol.models.classification\_tools module

Module containing functions for performing classification, via machine learning models or otherwise, related to Brain-Computer Interface applications.

`neurol.models.classification_tools.get_channels` (*signal*, *channels*, *device=None*)

Returns a signal with only the desired channels.

#### Parameters

- **signal** (*np.ndarray*) – a signal of shape [n\_samples, n\_channels]
- **channels** (*array*) – str names or int indices of the desired channels. returned in given order.
- **device** (*str*) – name of the device. Optional.

**Returns** numpy array of signal with shape [n\_channels, n\_desired\_channels]. Includes only the selected channels in the order given.

`neurol.models.classification_tools.softmax_predict` (*input\_*, *predictor*, *thresh=0.5*)

Consolidates a softmax prediction to a one-hot encoded prediction.

#### Parameters

- **input** – the input taken by the predictor
- **predictor** – function which returns a softmax prediction given an **input\_**
- **thresh** – the threshold for a positive prediction for a particular class.

`neurol.models.classification_tools.encode_ohe_prediction` (*prediction*)

Returns the index number of the positive class in a one-hot encoded prediction.

`neuro1.models.classification_tools.decode_prediction(prediction, decode_dict)`

Returns a more intelligible reading of the prediction based on the given `decode_dict`

`neuro1.models.classification_tools.threshold_clf(features, threshold, clf_consolidator='any')`

Classifies given features based on a given threshold.

#### Parameters

- **features** – an array of numerical features to classify
- **threshold** – threshold for classification. A single number, or an array corresponding to *features* for element-wise comparison.
- **clf\_consolidator** – method of consolidating element-wise comparisons with threshold into a single classification.  
 'any': positive class if any features passes the threshold  
 'all': positive class only if all features pass threshold  
 'sum': a count of the number of features which pass the threshold  
 function: a custom function which takes in an array of booleans and returns a consolidated classification

**Returns** classification for the given features. Return type *clf\_consolidator*.

### 1.1.2 neuro1.models.data\_exploration module

Module containing functions to study and analyze neural signals, especially to provide insights for building machine learning models to perform classification relevant to Brain-Computer Interface applications.

`neuro1.models.data_exploration.plot_signal(signal, sampling_rate, signal_type=None, ch_names=None, event_timestamps=None, **plt_kwargs)`

Plots signal.

#### Parameters

- **signal** – signal as array of shape [n\_samples, n\_channels].
- **sr** (*float*) – sampling rate in samples per second.
- **signal\_type** – (optional) gives a title for the y-axis.
- **ch\_names** – (optional) array of names for each channel (used for legend).
- **event\_timestamps** – (optional) 1-D array of times at which an event/stimulus occurred.
- **\*\*plt\_kwargs** – matplotlib keyword args

`neuro1.models.data_exploration.plot_grid(signals, num_signals=None, sampling_rate=1, cols=4, fig_size=(10, 6), sharey=True, sharex=True, random=True, fig_axes=None, show=True)`

Plot an (optionally random) set of signals [epochs] in a grid from a larger array of given signals.

#### Parameters

- **signals** – array of signals to plot from (num\_signals, num\_samples).
- **num\_signals** (*int*) – the number of signals to plot.
- **sampling\_rate** (*float*) – sampling rate of signals.
- **cols** (*int*) – the number of columns in the grid.



- **fig\_size** – tuple (x,y) of figure size in inches.
- **sharey** (*bool*) – whether to share scale on y-axis (see matplotlib).
- **sharex** (*bool*) – whether to share scale on x-axis (see matplotlib).
- **random** (*bool*) – whether to choose signals randomly or just use the first num\_signals.
- **fig\_axes** – optionally, an existing tuple of (fig,axes) to plot on (see matplotlib) rather creating new one.
- **show** (*bool*) – whether to show plot inline.

**Returns**

matplotlib figure and axes with sample of signals plotted in a grid

**Return type** fig, axes

neuro1.models.data\_exploration.**stim\_triggered\_average** (*signal, sampling\_rate, timestamps, duration\_before, duration\_after, plot=False*)

Inspired by the computational neuroscience concept of the spike-triggered average, this function computes the average signal characteristic around known events.

**Parameters**

- **signal** – signal as an array of shape [samples, channels].
- **sr** (*float*) – sampling rate of the signal.
- **timestamps** – array of floats containing the timestamps for each event.
- **duration\_before** – the duration to be considered before each event.
- **duration\_after** – the duration to be considered after each event.
- **plot** (*optional*) – whether or not to plot the stim\_triggered\_average.

**Returns**

average signal characteristic around event. relative\_time: relative time of each sample in stim\_triggered\_average  
with respect to event.

**Return type** stim\_triggered\_average

neuro1.models.data\_exploration.**plot\_PCA** (*epochs, sampling\_rate=1, n\_components=None, return\_PCA=False, PCA\_kwargs=None, plot\_grid\_kwargs=None*)

performs principal component analysis and plots principal components of epochs of a signal.

**Parameters**

- **epochs** – array of epochs (n\_epochs, n\_samples).
- **sr** (*float*) – sampling rate.
- **num\_components** (*int*) – number of components to use. If none is passed, all are used.
- **return\_PCA** (*bool*) – whether to return the independent components.
- **PCA\_kwargs** (*dict*) – dictionary containing kwargs for PCA function (see scikit-learn).
- **plot\_grid\_kwargs** (*dict*) – dictionary containing kwargs for plot\_grid function.

`neuro1.models.data_exploration.plot_ICA`(*epochs*, *sampling\_rate=1*, *n\_components=None*,  
*return\_ICA=False*, *FastICA\_kwargs=None*,  
*plot\_grid\_kwargs=None*)

performs independent component analysis and plots independent components of epochs of a signal.

#### Parameters

- **epochs** – array of epochs (*n\_epochs*, *n\_samples*).
- **sr** (*float*) – sampling rate.
- **num\_components** (*int*) – number of components to use. If none is passed, all are used.
- **return\_ICA** (*bool*) – whether to return the independent components.
- **FastICA\_kwargs** (*dict*) – dictionary containing kwargs for FastICA function (see scikit-learn).
- **plot\_grid\_kwargs** (*dict*) – dictionary containing kwargs for plot\_grid function.

### 1.1.3 neuro1.models.model\_tools module

Module for managing the models which come pre-packaged with neuro1. Includes functionality for importing and using the models.

`neuro1.models.model_tools.get_model`(*model\_name*)

gets the specified trained model.

**Parameters** *model\_name* (*str*) – name of model. See documentation for list of available models.

**Returns** trained model.

**Return type** model

`neuro1.models.model_tools.get_predictor`(*model\_name*)

gets the predictor for the specified model.

**Parameters** *model\_name* (*str*) – name of model. See documentation for list of available models.

**Returns** predictor of trained model.

**Return type** predictor

### 1.1.4 neuro1.models.preprocessing module

Module containing functions for the preparation of neural data for use with with BCI-related models.

`neuro1.models.preprocessing.epoch`(*signal*, *window\_size*, *inter\_window\_interval*)

Creates overlapping windows/epochs of EEG data from a single recording.

#### Parameters

- **signal** – array of timeseries EEG data of shape [*n\_samples*, *n\_channels*]
- **window\_size** (*int*) – desired size of each window in number of samples
- **inter\_window\_interval** (*int*) – interval between each window in number of samples (measured start to start)

**Returns** numpy array object with the epochs along its first dimension

`neuro1.models.preprocessing.labels_from_timestamps` (*timestamps*, *sampling\_rate*, *length*)

takes an array containing timestamps (as floats) and returns a labels array of size 'length' where each index corresponding to a timestamp via the 'samplingRate'.

#### Parameters

- **timestamps** – an array of floats containing the timestamps for each
- **event** (*units matching sampling\_rate*) –
- **sampling\_rate** (*float*) – the sampling rate of the EEG data.
- **length** (*int*) – the number of samples of the corresponding EEG recording.

**Returns** an integer array of size 'length' with a '1' at each time index where a corresponding timestamp exists, and a '0' otherwise.

`neuro1.models.preprocessing.label_epochs` (*labels*, *window\_size*, *inter\_window\_interval*, *label\_method*)

create labels for individual epochs of EEG data based on the label\_method.

#### Parameters

- **labels** – an integer array indicating a class for each sample measurement
- **window\_size** (*int*) – size of each window in number of samples (matching window\_size in epoched data)
- **inter\_window\_interval** (*int*) – interval between each window in number of samples (matching inter\_window\_interval in epoched data)
- **label\_method** (*str/func*) – method of consolidating labels contained in epoch into a single label.

'containment': whether a positive label occurs in the epoch, 'count': the count of positive labels in the epoch, 'mode': the most common label in the epoch func: func\_name(epoched\_labels) outputs label of epoched\_labels

**Returns** a numpy array with a label corresponding to each epoch

`neuro1.models.preprocessing.label_epochs_from_timestamps` (*timestamps*, *sampling\_rate*, *length*, *window\_size*, *inter\_window\_interval*, *label\_method='containment'*)

Directly creates labels for individual windows of EEG data from timestamps of events.

#### Parameters

- **timestamps** – an array of floats containing the timestamps for each event (units matching sampling\_rate).
- **sampling\_rate** (*float*) – sampling rate of the recording.
- **length** (*int*) – the number of samples of the corresponding EEG recording.
- **window\_size** (*int*) – size of each window in number of samples (matches window\_size in epoched data)
- **inter\_window\_interval** (*int*) – interval between each window in number of samples (matches inter\_window\_interval in epoched data)
- **label\_method** (*str/func*) – method of consolidating labels contained in epoch into a single label.

'containment': whether a positive label occurs in the epoch, 'count': the count of positive labels in the epoch, 'mode': the most common label in the epoch func: func\_name(epoched\_labels) outputs label of epoched\_labels

**Returns** an array with a label correponding to each window

neuro1.models.preprocessing.**epoch\_and\_label**(data, sampling\_rate, timestamps, window\_size, inter\_window\_interval, label\_method='containment')

Epochs a signal (single EEG recording) and labels each epoch using timestamps of events and a chosen labelling method.

#### Parameters

- **data** – array of timeseries EEG data of shape [n\_samples, n\_channels]
- **timestamps** – an array of floats containing the timestamps for each event in units of time.
- **sampling\_rate** (*float*) – the sampling rate of the EEG data.
- **window\_size** (*float*) – desired size of each window in units of time.
- **inter\_window\_interval** (*float*) – interval between each window in units of time (measured start to start)
- **label\_method** (*str/func*) – method of consolidating labels contained in epoch into a single label.

'containment': whether a positive label occurs in the epoch, 'count': the count of positive labels in the epoch, 'mode': the most common label in the epoch func: func\_name(epoched\_labels) outputs label of epoched\_labels

**Returns** array of epochs with shape [n\_epochs, n\_channels] labels: array of labels corresponding to each epoch of shape [n\_epochs,]

**Return type** epochs

neuro1.models.preprocessing.**compute\_signal\_std**(signal, corrupt\_intervals=None, sampling\_rate=1)

Computes and returns the standard deviation of a signal channel-wise while avoiding corrupt intervals

#### Parameters

- **signal** – signal of shape [n\_samples, n\_channels]
- **corrupt\_intervals** – an array of 2-tuples indicating the start and end time of the corrupt interval (units of time)
- **sampling\_rate** – the sampling rate in units of samples/unit of time

**Returns** standard deviation of signal channel-wise of shape [1, n\_channels]

neuro1.models.preprocessing.**split\_corrupt\_signal**(signal, corrupt\_intervals, sampling\_rate=1)

Splits a signal with corrupt intervals and returns array of signals with the corrupt intervals filtered out. This is useful for treating each non\_corrupt segment as a seperate signal to ensure continuity within a single signal.

#### Parameters

- **signal** – signal of shape [n\_samples, n\_channels]
- **corrupt\_intervals** – an array of 2-tuples indicating the start and end time of the corrupt interval (units of time)
- **sampling\_rate** – the sampling rate in units of samples/unit of time

**Returns** array of non\_corrupt signals of shape [n\_signal, n\_samples, n\_channels]

neuroI.models.preprocessing.**epoch\_band\_features**(*epoch\_*, *sampling\_rate*, *bands*='all',  
return\_dict=True)

Computes power features of EEG frequency bands over the epoch channel-wise.

#### Parameters

- **epoch** – a single epoch of shape [n\_samples, n\_channels]
- **sampling\_rate** – the sampling rate of the signal in units of samples/sec
- **bands** – the requested frequency bands to get power features for. 'all': all of ['theta', 'alpha\_low', 'alpha\_high', 'beta', 'gamma'] otherwise an array of strings of the desired bands.
- **return\_dict** (*bool*) – returns *band\_features* in the form of a dictionary if True, else returns as numpy array in order of *bands*

**Returns** a dictionary of arrays of shape [1, n\_channels] containing the power features over each frequency band per channel.

## 1.2 neuroI.BCI module

Module implementing a general Brain-Computer Interface which manages and incoming stream of neural data and responds to it in real-time.

**class** neuroI.BCI.**generic\_BCI**(*classifier*, *transformer*=None, *action*=<built-in function print>,  
calibrator=None)

Bases: object

Implements a generic Brain-Computer Interface.

Internally manages a buffer of the signal given in *stream* and continuously performs classification on appropriately transformed real-time neural data. At each classification, a corresponding *action* is performed.

#### Variables

- **classifier** (*function*) – a function which performs classification on the most recent data (transformed as needed). returns classification.
- **transformer** (*function*) – function which takes in the most recent data (*buffer*) and returns the transformed input the classifier expects.
- **action** (*function*) – a function which takes in the classification, and performs some action.
- **calibrator** (*function*) – a function which is run on startup to perform calibration using *stream*; returns *calibration\_info* which is used by *classifier* and *transformer*.
- **calibration\_info** – the result of the calibrator, if applicable.
- **buffer\_length** (*int*) – the length of the *buffer*; specifies the number of samples of the signal to keep for classification.
- **brain\_state** – the most recent brain state classification.

**\_\_init\_\_**(*classifier*, *transformer*=None, *action*=<built-in function print>, *calibrator*=None)

Initialize a generic BCI object.

See class documentation for information about the class itself.

#### Parameters

- **classifier** (*function*) – a function which performs classification on the most recent data (transformed as needed). returns class.
- **transformer** (*function*) – function which takes in the most recent data (*buffer*) and returns the transformed input the classifier expects.
- **action** (*function*) – a function which takes in the classification, and performs some action.
- **calibrator** (*function*) – a function which is run on startup to perform calibration using *stream*; returns *calibration\_info* which is used by *classifier* and *transformer*.

**calibrate** (*stream*)

runs the *calibrator*.

return value of *calibrator* is stored in the object's *calibration\_info* which the *transformer* and *classifier* can use at run-time of BCI.

**Parameters** **stream** (*neuro1.streams object*) – neuro1 stream for brain data.

**run** (*stream*)

Runs the defined Brain-Computer Interface.

Internally manages a buffer of the signal given in *stream* and continuously performs classification on appropriately transformed real-time neural data. At each classification, a corresponding *action* is performed.

**Parameters** **stream** (*neuro1.streams object*) – neuro1 stream for brain data.

**test\_update\_rate** (*stream, test\_length=10, perform\_action=True*)

Returns the rate at which the BCI is able to make a classification and perform its action.

**Parameters**

- **stream** (*neuro1.streams object*) – neuro1 stream for brain data.
- **test\_length** (*float*) – how long to run the test for in seconds.
- **perform\_action** (*bool*) – whether to perform the action or skip it.

**class** neuro1.BCI.fsm\_BCI (*classifier, transformer=None, action=<built-in function print>, calibrator=None*)

Bases: neuro1.BCI.generic\_BCI

Implements a Finite-State-Machine-inspired Brain-Computer Interface.

Classification of brain-state is not only dependent on the transformed real-time brain signal, but also the previous brain state.

Internally manages a buffer of the signal given in *stream* and continuously performs classification on appropriately transformed real-time neural data. At each classification, a corresponding *action* is performed.

**Variables**

- **classifier** (*function*) – a function which performs classification on the most recent data (transformed as needed). returns classification.
- **transformer** (*function*) – function which takes in the most recent data (*buffer*) and returns the transformed input the classifier expects.
- **action** (*function*) – a function which takes in the classification, and performs some action.
- **calibrator** (*function*) – a function which is run on startup to perform calibration using *stream*; returns *calibration\_info* which is used by *classifier* and *transformer*.
- **calibration\_info** – the result of the calibrator, if applicable.

- **buffer\_length** (*int*) – the length of the *buffer*; specifies the number of samples of the signal to keep for classification.
- **brain\_state** – the most recent brain state classification.

**class** neuro1.BCI.**retentive\_BCI** (*classifier*, *transformer*=None, *action*=<built-in function print>, *calibrator*=None, *memory\_length*=10)

Bases: neuro1.BCI.generic\_BCI

Implements a Brain-Computer Interface with memory of past brain states.

Classification of brain-state is not only dependent on the transformed real-time brain signal, but also the finite list of previous brain states.

Internally manages a buffer of the signal given in *stream* and continuously performs classification on appropriately transformed real-time neural data. At each classification, a corresponding *action* is performed.

#### Variables

- **classifier** (*function*) – a function which performs classification on the most recent data (transformed as needed). returns classification.
- **transformer** (*function*) – function which takes in the most recent data (*buffer*) and returns the transformed input the classifier expects.
- **action** (*function*) – a function which takes in the classification, and performs some action.
- **calibrator** (*function*) – a function which is run on startup to perform calibration using *stream*; returns *calibration\_info* which is used by *classifier* and *transformer*.
- **calibration\_info** – the result of the calibrator, if applicable.
- **brain\_state** – the most recent brain state classification.
- **memory\_length** (*int*) – number of brain states into the past to remember.
- **past\_states** – a list of the past classifications of brain states. used in next classification. length is *memory\_length*.

**\_\_init\_\_** (*classifier*, *transformer*=None, *action*=<built-in function print>, *calibrator*=None, *memory\_length*=10)

Initialize a retentive BCI object.

See class documentation for information about the class itself.

#### Parameters

- **classifier** (*function*) – a function which performs classification on the most recent data (transformed as needed). returns class.
- **transformer** (*function*) – function which takes in the most recent data (*buffer*) and returns the transformed input the classifier expects.
- **action** (*function*) – a function which takes in the classification, and performs some action.
- **calibrator** (*function*) – a function which is run on startup to perform calibration using *stream*; returns *calibration\_info* which is used by *classifier* and *transformer*.
- **memory\_length** (*int*) – number of brain states to remember into past.

**class** neuro1.BCI.**automl\_BCI** (*model*, *epoch\_len*, *n\_states*, *transformer*=None, *action*=<built-in function print>)

Bases: neuro1.BCI.generic\_BCI

Implements a Brain-Computer Interface which builds its own classifier by training a machine learning model in the calibration stage.

At calibration, data is recorded for some number of brain-states then a machine-learning classifier is trained on the transformed data.

Internally manages a buffer of the signal given in *stream* and continuously performs classification on appropriately transformed real-time neural data. At each classification, a corresponding *action* is performed.

#### Variables

- **model** – a model object which has `fit(X, y)` and `predict(X)` methods.
- **classifier** (*function*) – the model’s predictor after training. accepts transformed data and returns classification.
- **transformer** (*function*) – function which takes in the most recent data (*buffer*) and returns the transformed input the classifier expects.
- **action** (*function*) – a function which takes in the classification, and performs some action.
- **brain\_state** – the most recent brain state classification.

**\_\_init\_\_** (*model, epoch\_len, n\_states, transformer=None, action=<built-in function print>*)  
Initialize an autoML BCI object.

See class documentation for information about the class itself.

#### Parameters

- **model** – a model object which has `fit(X, y)` and `predict(X)` methods.
- **epoch\_len** (*int*) – the length of the epochs (in # of samples) used in training and prediction by the model.
- **n\_states** (*int*) – the number of brain states being classified.
- **transformer** (*function, optional*) – function which takes in the most recent data (*buffer*) and returns the transformed input the classifier expects. Defaults to None.
- **action** (*function, optional*) – a function which takes in the classification, and performs some action. Defaults to print.

**build\_model** (*stream, recording\_length*)  
records brain signal

#### Parameters

- **stream** (*neuro1.streams object*) – neuro1 stream for brain data.
- **recording\_length** (*float*) – length in seconds for the recording of each brain state to be used for training the model.

**run** (*stream*)  
Runs the defined Brain-Computer Interface.

Internally manages a buffer of the signal given in *stream* and continuously performs classification on appropriately transformed real-time neural data. At each classification, a corresponding *action* is performed.

**Parameters** **stream** (*neuro1.streams object*) – neuro1 stream for brain data.



## 1.3 neuro1.BCI\_tools module

Module including utility functions for creating *classifier*'s, *transformer*'s, and *calibrator*'s for use in the *BCI* module.

`neuro1.BCI_tools.ensemble_transform` (*signal*, *epoch\_len=None*, *channels=None*, *device=None*,  
*transformers=None*, *filter\_=False*, *sampling\_rate=None*,  
*filter\_kwargs=None*)

Ensemble transform function. Takes in buffer as input. Extracts the appropriate channels and samples, performs filtering, and transforms.

### Parameters

- **signal** (*np.ndarray*) – signal of shape: [n\_samples, n\_channels]
- **epoch\_len** (*int*) – length of epoch expected by classifier (# of samples). optional.
- **channels** (*list of str or int*) – list of channels expected by classifier. See `get_channels`. optional.
- **device** (*str*) – device name. used to get channels and *sampling\_rate*.
- **filter** (*boolean*) – whether to perform filtering
- **filter\_kwargs** (*dict*) – dictionary of kwargs passed to filtering function. See `biosppy.signals.tools.filter_signal`. by default, an order 8 bandpass butter filter is performed between 2Hz and 40Hz.

`neuro1.BCI_tools.filter_signal` (*signal*, *sampling\_rate*, *ftype='butter'*, *band='bandpass'*, *frequency=(2, 40)*, *order=8*, *\*\*filter\_kwargs*)

applies frequency-based filters to a given signal.

### Parameters

- **signal** (*np.ndarray*) – signal of shape [n\_samples, n\_channels]
- **sampling\_rate** (*float*) – sampling rate of signal.
- **ftype** (*str, optional*) – type of filter. one of 'FIR', 'butter', 'chebby1', 'chebby2', 'ellip', or 'bessel'. Defaults to 'butter'.
- **band** (*str, optional*) – band type. one of 'lowpass', 'highpass', 'bandpass', or 'band-stop'. Defaults to 'bandpass'.
- **frequency** (*float or tuple of floats, optional*) – cutoff frequencies. single if 'lowpass'/'highpass', tuple if 'bandpass'/'bandstop'. Defaults to (2,40).
- **order** (*int, optional*) – order of filter. Defaults to 8.
- **\*\*filter\_kwargs** – keyword args for `biosppy.signals.tools.filter_signal`

**Returns** filtered signal

**Return type** [np.ndarray]

`neuro1.BCI_tools.band_power_calibrator` (*stream*, *channels*, *device*, *bands*, *percentile=50*,  
*recording\_length=10*, *epoch\_len=1*, *inter\_window\_interval=0.2*)

Calibrator for *generic\_BCI.BCI* which computes a given *percentile* for the power of each frequency band across epochs channel-wise. Useful for calibrating a concentration-based BCI.

### Parameters

- **stream** (*neuro1.streams object*) – neuro1 stream for brain data.
- **channels** – array of strings with the names of channels to use.

- **device** (*str*) – device name for use by *classification\_tools*
- **bands** – the frequency bands to get power features for. ‘all’: all of [‘theta’, ‘alpha\_low’, ‘alpha\_high’, ‘beta’, ‘gamma’] otherwise an array of strings of the desired bands.
- **percentile** – the percentile of power distribution across epochs to return for each band.
- **recording\_length** (*float*) – length of recording to use for calibration in seconds.
- **epoch\_len** (*float*) – the length of each epoch in seconds.
- **inter\_window\_interval** (*float*) – interval between each window/epoch in seconds (measured start to start)

**Returns** array of shape [n\_bands, n\_channels] of the *percentile* of the power of each band

**Return type** clb\_info

`neuroI.BCI_tools.band_power_transformer` (*signal, sampling\_rate, bands*)

Transformer for *generic\_BCI.BCI* which chooses channels, epochs, and gets power features on some choice of bands.

**Parameters**

- **signal** (*np.ndarray*) – most recent stream data. shape: [n\_samples, n\_channels]
- **sampling\_rate** (*float*) – sampling\_rate of signal.
- **bands** – the frequency bands to get power features for. ‘all’: all of [‘theta’, ‘alpha\_low’, ‘alpha\_high’, ‘beta’, ‘gamma’] otherwise a list of strings of the desired bands.

**Returns** array of shape [n\_bands, n\_channels] of the channel-wise power of each band over the epoch.

**Return type** transformed\_signal

## 1.4 neuroI.plot module

Module for plotting stream of neural data. Includes time domain, fourrier transform, and spectrogram live plots.

`neuroI.plot.plot` (*stream, channels=None, w\_size=(1920, 1080)*)

plots data stream. one row per channel.

**Parameters**

- **stream** (*neuroI.streams object*) – neuroI stream for a data source.
- **channels** – channels to plot. list/tuple of channel indices, or dict with indices as keys and names as values. Defaults to None (plots all channels w/o names).
- **w\_size** (*tuple, optional*) – initial size of window in pixels. Defaults to (1920, 1080).

`neuroI.plot.plot_fft` (*stream, channels=None, w\_size=(1920, 1080)*)

plots fourrier transform of data stream from inlet. one row per channel.

**Parameters**

- **stream** (*neuroI.streams object*) – neuroI stream for a data source.
- **channels** – channels to plot. list/tuple of channel indices, or dict with indices as keys and names as values. Defaults to None (plots all channels w/o names).

- **w\_size** (*tuple, optional*) – initial size of window in pixels. Defaults to (1920, 1080).

`neuro1.plot.plot_spectrogram(stream, channels=None, w_size=(1920, 1080))`  
 plots spectrogram of data stream from inlet. one row per channel.

#### Parameters

- **stream** (*neuro1.streams object*) – neuro1 stream for a data source.
- **channels** – channels to plot. list/tuple of channel indices, or dict with indices as keys and names as values. Defaults to None (plots all channels w/o names).
- **w\_size** (*tuple, optional*) – initial size of window in pixels. Defaults to (1920, 1080).

## 1.5 neuro1.connect\_device module

Module containing functions for quickly connecting to BCI-related streaming devices.

`neuro1.connect_device.connect_muse()`  
 connects to any available muse headset. returns `ble2lsl.ble2lsl.Streamer` object.

`neuro1.connect_device.get_lsl_EEG_inlets()`  
 resolves all EEG lsl streams and returns their inlets in an array.

## 1.6 neuro1.streams

module for handling streams of data from different sources

**class** `neuro1.streams.lsl_stream` (*pylsl\_inlet, buffer\_length=2048*)  
 Bases: `object`

A generalized stream object for an lsl data source.

Manages a buffer of data and makes it available. Used by `neuro1.BCI` and `neuro1.plot`.

**\_\_init\_\_** (*pylsl\_inlet, buffer\_length=2048*)  
 initialize an `lsl_stream` object.

#### Parameters

- **pylsl\_inlet** (*pylsl.pylsl.StreamInlet*) – inlet of connected lsl device
- **buffer\_length** (*int, optional*) – length of data buffer. Defaults to 2048.

**get\_data** (*max\_samples=2048*)  
 gets latest data.

**record\_data** (*duration*)  
 records from stream for some duration of time.

**Parameters** **duration** (*float*) – length of recording in seconds.

**update\_buffer** ()  
 updates buffer with most recent available data.

**Returns** True if new data available, False if not.

**Return type** [bool]

**close()**  
closes the pylsl inlet stream

## CHAPTER 2

---

### Installation

---

neuro1 can be easily installed using pip:

```
$ pip install neuro1
```



## CHAPTER 3

---

### Indices

---

- `genindex`
- `modindex`
- `search`





## Symbols

[\\_\\_init\\_\\_\(\) \(neuroI.BCI.automI\\_BCI method\), 12](#)  
[\\_\\_init\\_\\_\(\) \(neuroI.BCI.generic\\_BCI method\), 9](#)  
[\\_\\_init\\_\\_\(\) \(neuroI.BCI.retentive\\_BCI method\), 11](#)  
[\\_\\_init\\_\\_\(\) \(neuroI.streams.lsl\\_stream method\), 15](#)

## A

[automI\\_BCI \(class in neuroI.BCI\), 11](#)

## B

[band\\_power\\_calibrator\(\) \(in module neuroI.BCI\\_tools\), 13](#)  
[band\\_power\\_transformer\(\) \(in module neuroI.BCI\\_tools\), 14](#)  
[build\\_model\(\) \(neuroI.BCI.automI\\_BCI method\), 12](#)

## C

[calibrate\(\) \(neuroI.BCI.generic\\_BCI method\), 10](#)  
[close\(\) \(neuroI.streams.lsl\\_stream method\), 15](#)  
[compute\\_signal\\_std\(\) \(in module neuroI.models.preprocessing\), 8](#)  
[connect\\_muse\(\) \(in module neuroI.connect\\_device\), 15](#)

## D

[decode\\_prediction\(\) \(in module neuroI.models.classification\\_tools\), 3](#)

## E

[encode\\_ohe\\_prediction\(\) \(in module neuroI.models.classification\\_tools\), 3](#)  
[ensemble\\_transform\(\) \(in module neuroI.BCI\\_tools\), 13](#)  
[epoch\(\) \(in module neuroI.models.preprocessing\), 6](#)  
[epoch\\_and\\_label\(\) \(in module neuroI.models.preprocessing\), 8](#)  
[epoch\\_band\\_features\(\) \(in module neuroI.models.preprocessing\), 9](#)

## F

[filter\\_signal\(\) \(in module neuroI.BCI\\_tools\), 13](#)  
[fsm\\_BCI \(class in neuroI.BCI\), 10](#)

## G

[generic\\_BCI \(class in neuroI.BCI\), 9](#)  
[get\\_channels\(\) \(in module neuroI.models.classification\\_tools\), 3](#)  
[get\\_data\(\) \(neuroI.streams.lsl\\_stream method\), 15](#)  
[get\\_lsl\\_EEG\\_inlets\(\) \(in module neuroI.connect\\_device\), 15](#)  
[get\\_model\(\) \(in module neuroI.models.model\\_tools\), 6](#)  
[get\\_predictor\(\) \(in module neuroI.models.model\\_tools\), 6](#)

## L

[label\\_epochs\(\) \(in module neuroI.models.preprocessing\), 7](#)  
[label\\_epochs\\_from\\_timestamps\(\) \(in module neuroI.models.preprocessing\), 7](#)  
[labels\\_from\\_timestamps\(\) \(in module neuroI.models.preprocessing\), 6](#)  
[lsl\\_stream \(class in neuroI.streams\), 15](#)

## N

[neuroI.BCI \(module\), 9](#)  
[neuroI.BCI\\_tools \(module\), 13](#)  
[neuroI.connect\\_device \(module\), 15](#)  
[neuroI.models.classification\\_tools \(module\), 3](#)  
[neuroI.models.data\\_exploration \(module\), 4](#)  
[neuroI.models.model\\_tools \(module\), 6](#)  
[neuroI.models.preprocessing \(module\), 6](#)  
[neuroI.plot \(module\), 14](#)  
[neuroI.streams \(module\), 15](#)

## P

[plot\(\) \(in module neuroI.plot\), 14](#)

`plot_fft()` (*in module neurol.plot*), 14  
`plot_grid()` (*in module neurol.models.data\_exploration*), 4  
`plot_ICA()` (*in module neurol.models.data\_exploration*), 5  
`plot_PCA()` (*in module neurol.models.data\_exploration*), 5  
`plot_signal()` (*in module neurol.models.data\_exploration*), 4  
`plot_spectrogram()` (*in module neurol.plot*), 15

## R

`record_data()` (*neurol.streams.lsl\_stream method*), 15  
`retentive_BCI` (*class in neurol.BCI*), 11  
`run()` (*neurol.BCI.automl\_BCI method*), 12  
`run()` (*neurol.BCI.generic\_BCI method*), 10

## S

`softmax_predict()` (*in module neurol.models.classification\_tools*), 3  
`split_corrupt_signal()` (*in module neurol.models.preprocessing*), 8  
`stim_triggered_average()` (*in module neurol.models.data\_exploration*), 5

## T

`test_update_rate()` (*neurol.BCI.generic\_BCI method*), 10  
`threshold_clf()` (*in module neurol.models.classification\_tools*), 4

## U

`update_buffer()` (*neurol.streams.lsl\_stream method*), 15